

# Installing `p1z`

---

Basically, `p1z.exe` is a completely self-contained executable. Just copy it onto your machine and double-click... This will *not* install any software permanently -- it will simply run. That means you don't need to bother keeping it on your machine, just get another copy next time you need to run a pipeline on your PC!

The big exception is when running `dewey`, which is the library search command. For this you will need to have three libraries on your machine: NIST msms, Metlin and our very own Lipid\_MBX.

Simply place the folders (found in `Official_Spectral_Libraries` ) named:

- `nist_msms`
- `METLIN_EXPERIMENTAL`
- `LipidBlast_MBX`
- `DECOY_nist_msms`
- `DECOY_METLIN_EXPERIMENTAL`
- `DECOY_LipidBlast_MBX`

Into the following directory: `C:\NIST14\MSSEARCH`